Overview ty-outlier-explorer - Quality Control Details

This page provides a detailed report of the quality control documentation for safety-outlier-explorer. The source code for this page is available on github.

Technical Specs

Technical documentation including, functional specifications, regressions tests and risk assessments. see original

Functional Specifications

Filter by Measure

This drop-down menu is used to filter the lab variable that is displayed in the chart. The default measure is set to the first variable in the data set, but the default value can be adjusted to any measure upon the project team's request. Users can select any pre-defined measure in this filter. Selecting a measure will display the data points associated with that measure and will display the visits that the measure is collected at underneath the chart.

Filter by pre-selected characteristics

The chart can be configured to have any number of data filters. Each specified filter has a drop-down menu that is used to filter the visit data that is displayed in the chart.

Total Participants Shown

Displayed under the header, this statement provides the total number of participants in the data set in integers as well as the percentage of participants shown to one decimal place. In the default view, all participant data is included in the chart, or 100%, unless the project team requests otherwise. When filters are used, the data shown in the histogram is manipulated and the number of participant records shown at the top of the page may increase or decrease, depending on the variable selected.

Change the X-axis

The chart contains a toggle that controls the display of the x-axis. Users can view the x-axis data by study day, visit, or visit number. The x-axis displays the data point at which the measure data is collected and manipulating it changes the scale at which the measure data is plotted. Study day, visit, or visit number can be chosen for a given measure, but the chart default is to show the x-axis by study day.

Set Y-axis Limits

Users are able to adjust the lower limit of the y-axis by manipulating the input within the Lower Limit box. Users are able to adjust the upper limit of the y-axis by manipulating the input within the Upper Limit box. Each limit box has up and down arrows to adjust the value of the box by 1. Users can select an arrow to change the input or delete the current value in the box and type in a new value. When typing in a value, users must click outside of the box to apply the changes. When a user manipulates a y-axis limit, the y-axis on the chart below the controls will reflect the change.

Reset Y-axis limits

Users can change the upper and lower y-axis limits back to the default values for a measure by selecting the "Reset Limits" button. Pressing the button will erase the current values in the upper and lower limit boxes and replace them with defaults. The chart below will also reflect the changes and will show the

default axes.

Normal Range Controls

Users can opt to change the method for the normal ranges to the following options: None, LLN-ULN, Standard Deviation, and Quantiles. The default for the normal range controls is set to "LLN-ULN", which is the upper and lower normal range specified in the data. Users may select "None", which will remove the green shading around the normal ranges and will remove the inlier count to the left of the charts. Users may also select "Standard Deviation", which will prompt the user to enter the number of standard deviations they would like to set as the normal range. Altering the default value will change the upper and lower normal range bands in the chart and will change the inlier count to the left of the chart. Users may also select "Quantiles", which prompts the user to enter an upper and lower quantile to set as the normal ranges. Altering the upper quantile will change the upper normal range band in the chart and altering the lower quantile will change the lower normal range band in the chart. Changing the normal ranges may change the number of inliers and the inlier count to the left of the chart will update based on your selections.

View Visits Without Data

Users can select the checkbox to view visits that do not contain data. Checking the box will show a timepoint for the visit. Unchecking the checkbox will hide the timepoint for visits without data; visits without data are hidden on default.

View Unscheduled Visits

Users can select the checkbox to view data that is not associated with a scheduled visit. Checking the box will show a timepoint for the unscheduled visit and its associated data. Unselecting the checkbox will hide the unscheduled visit timepoint and its data; unscheduled visits are hidden on default.

View Small Multiples

Clicking a data point will open a listing of small multiples below the chart that contains a detailed listing of all of the participant's records. This linked table contains small charts for each measure captured in the data set as well as the participant's subject ID, age, sex, and race. When the listing appears, the data points selected will appear in a continuous gray line inside of the chart. When a user clicks outside of the selected data point, the gray line disappears along with the detailed listing.

Clicking a data point will perform 2 actions:

- Display the raw data records, including the subject identifier, age, race, sex, and small data renderings of each measure captured in the data set
- Create a gray line to connect the selected data points

Regression Tests

- Confirm that the number and percentage of participants is captured above the chart
- Confirm that selecting the Measure filter changes the distribution of observations
- Confirm that the x-axis filter changes the x-axis to the selected variable
- Confirm that entering a value into the Y-axis Lower Limit box and clicking outside of the box updates the lower limit of the y-axis
- Confirm that entering a value into the Upper Limit box and clicking outside of the box updates the upper limit of the y-axis
- Confirm that selecting the y-axis reset button returns the chart to the default values for that measure
- Confirm that checking the "Visits without data" checkbox shows a timepoint for visits with no associated data points
- Confirm that unchecking the "Visits without data" checkbox hides the timepoint for visits with no associated data points
- Confirm that checking the "Unscheduled visits" checkbox shows a timepont for unscheduled visits and their associated data
- Confirm that unchecking the "Unscheduled visits" checkbox hides the timepoint for unscheduled visits and their associated data

- Confirm that hovering on a data point displays tooltips for the data point with Subject, value, and visit
- Confirm that date can be added to the standard tooltip with these settings: {tooltipcols: [{label:"Date", valuecol:"DT"}]}
- Confirm that clicking on a data point applies a dark highlight to all of the participant's data points
- Confirm that clicking on a data point opens a detailed listing of a participant's data beneath the chart
- Confirm that when clicking a point that overlaps with another point.
- Confirm that all marks (points, lines) associated with that participant should highlight.
- Confirm that a note should appear beneath the chart detailing which participants have results overlapping with the click mark.
- Confirm that hovering over any ID in the note should highlight that participant in the chart above.
- Confirm that clicking any ID in the note should highlight that participant in the chart above and display that participant's measures below.
- Confirm that with an overlapping point highlighted, changing any control should remove the note but maintain the highlighting of the selected participant.
- Confirm that the selected participant is displayed in a dropdown with the current participant selected
- Confirm that the listing contains subject-level demographic data (Age, Sex, and Race appear by default, can add other Participant Details if wanted)
- Confirm that the listing contains subject-level information about measures and and visits
- Confirm that you can change the selected participant via the dropdown and the new subject should be highlighted in the main chart and the small multiples should update to reflect that participant's data.
- Confirm that a dropdown labeled Normal range method appears in the controls section with four options: None, LLN-ULN, Standard Deviation, and Quantiles.
- Confirm that selecting Standard Deviation displays a text input labeled Number of standard deviations and defaults to 1.96.
- Confirm that selecting Quantiles displays two text inputs, one for the lower quantile and one for the upper quantile, labeled Lower quantile and Upper

- quantile respectively, that default to .05 and .95.
- Confirm that the x-axis tick labels change rotation angles depending on the x-axis control variable.
- Confirm that a unit variable is appended to each measure in the Measure dropdown.
- Confirm that there is a specified unit variable in the form on the left panel of CAT.
- Confirm that there are values displayed next to each measure under the specified unit variable in the data set.
- Confirm that the chart lines, when the scrolled on, darken significantly.
 When the user removes his/her cursor from the chart line, it should return to its default state.
- Confirm that when clicking the arrows in the y-axis limit inputs, the limits should increment by a reasonable factor of the y-axis range, or roughly 1/15 of the range. For example, if the range is 150 the increment should be 10 and if the range is 1.5 the increment should be .1.
- Confirm that the minimum and maximum points are not clipped by the edge of the chart canvas and can be clicked.
- Confirm that hovering over a point that overlaps with one or more other points yields an additional line in the tooltip: Number of overlapping point(s) = n.
- Confirm that hovering over a point that does not overlap with any other points does not yield an additional line in the tooltip.

Warnings and Notes

- Confirm that the right side of the webpage between the controls and the chart a note should appear that details how many records with missing and/or non-numeric results were removed.
- Confirm that the test dataset already contains missing records. To test non-numeric records edit ADBDS.csv by adding non-numeric values to the STRESN column.
- Confirm that a warning appears in the console when an invalid initially displayed measure is entered, i.e. one that doesn't appear in the data. (Fill in some dummy text in the "Initially Displayed Measure" field in Settings)
- Verify that after entering the dummy text in the "Initially Displayed

- Measure" field and reloading the chart, that the Measure dropdown defaults to an existent measure
- Confirm that when entering a valid initially displayed measure but without units, (e.g. Temperature) that a warning appears in the console noting that the specified measure does not include units
- After entering the measure and reloading the chart, verify the Measure dropdown defaults to the specified measure with units

Highlighting

- Confirm that adding the following block to the Settings input will add custom marks to the chart: { "custom_marks": [{ "per": ["USUBJID", "VISIT", "TEST", "STRESN"], "type": "circle", "attributes": { "fill": "red", "fill-opacity": 1, "stroke": "black", "stroke-opacity": 1 }, "radius": 4, "values": { "DY": ["56", "112", "168", "224", "280", "336"] } }] }
- Verify that custom marks as well as default marks highlight as expected
- · Verify that custom marks as well as default marks can be clicked
- Verify that custom marks maintain their visual characteristics, e.g. fill, stroke, etc.

Data Guidelines

Data Specifications for the charts. see original

The Safety Outlier Explorer is initialized with JSON data files that match the format created by d3.csv(). The chart expects an ADaM-esque data structure by default but can be customized to work with any dataset that contains one row per participant per time point per measure and includes the required columns specified below.

Default Chart

The following columns are expected in a data file used to create an instance of the safety outlier explorer. A sample chart using this configuration is found here.

Data Specification

Settings Variable	Description	Default	Required?
id_col	Unique Subject Identifier	USUBJID	Υ
time <i>cols[0].value</i> col	Visit Name	VISIT	Υ
time <i>cols[0].order</i> col	Visit Name (N)	VISITNUM	
time <i>cols[1].value</i> col	Study Day of Visit	DY	
measure_col	Name of Measure or Test	TEST	Υ
value_col	Result or Finding	STRESN	Υ
unit_col	Units	STRESU	
normal <i>col</i> low	Normal Range Lower Limit	STRNLO	
normal <i>col</i> high	Normal Range Upper Limit	STRNHI	
details[0].value_col	Age	AGE	
details[1].value_col	Sex	SEX	

Settings

details[2].value_col Race

No customization is needed to initialize a chart when the data object contains all default variables above. Just pass an empty object ({}) for settings and initialize the chart like so: safetyOutlierExplorer('body', {}).init(data).

RACE

Customized Chart with Queries

The following data specs and chart settings can be used to create a chart with several custom filters, including Sex, Race and Query Status. A sample chart using this configuration is found here.

Data Specification

Settings Variable	Description	Default	Required?
id_col	Unique Subject Identifier	USUBJID	Υ

time <i>cols[0].value</i> col	Visit Name	VISIT	Υ
time <i>cols[0].order</i> col	Visit Name (N)	VISITNUM	
time <i>cols[1].value</i> col	Study Day of Visit	DY	
measure_col	Name of Measure or Test	TEST	Υ
value_col	Result or Finding	STRESN	Υ
unit_col	Units	STRESU	
normal <i>col</i> low	Normal Range Lower Limit	STRNLO	
normal <i>col</i> high	Normal Range Upper Limit	STRNHI	
details[0].value_col	Age	AGE	
details[1].value_col	Sex	SEX	
details[2].value_col	Race	RACE	
filter	Site ID	SITEID	Υ
filter	Treatment Arm	ARM	Υ
filter	Sex	SEX	Υ
filter	Race	RACE	Υ
filter	Query Flag ("Y" for yes)	QUERYFL	Υ
custom_marks[0].tooltip	Query Details	QUERYDETAILS	Y

Settings Object

```
value_col: 'DY',
        type: 'linear',
        label: 'Study Day',
        rotate_tick_labels: false,
        vertical_space: 0
    }
],
filters: [
    {value_col: 'ARM', label: 'Arm'},
    {value_col: 'SITEID', label: 'Site ID'},
    {value_col: 'SEX', label: 'Sex'},
    {value_col: 'RACE', label: 'Race'},
    {value_col: 'QUERYFL', label: 'Open Query?'}
],
details: Γ
    {value_col: 'ARM', label: 'Arm'},
    {value_col: 'AGE', label: 'Age'},
    {value_col: 'SEX', label: 'Sex'},
    {value_col: 'RACE', label: 'Race'}
],
custom_marks: [
    {
        per: ['USUBJID', 'VISIT', 'STRESN'],
        type: 'circle',
        attributes: {
            'fill-opacity': 1,
            'stroke': 'red',
            'fill': 'red'
        },
        values: {
            QUERYFL: ['Y']
        },
        tooltip: 'Query Details: [QUERYDETAILS]'
    }
],
multiples_sizing: {
```

```
width: 250,
    height: 75
}
```

API

Technical specifications for API. see original

safetyOutlierExplorer(element, settings)

a factory to create a custom Webcharts chart object

returns: chart

Param Type Description

element string CSS selector identifying the element in which to create the chart

settings object specifying options for how the chart is to settings object appear and behave. Options defined here overwrite default values; see Configuration

Events

participantsSelected

The custom participantsSelected event is dispatched to the overall chart wrapper (chart.wrap) whenever the details for a given participant are viewed (or cleared) by clicking on a line or point in the chart. The event has a custom data property holding an array with the selected ID when a new participant is selected ($\lceil "123-456-7" \rceil$) or an empty array when participant details are cleared.

Chart Configuration

Technical specifications for chart configuration. see original

The most straightforward way to customize the Safety Outlier Explorer is by using a configuration object whose properties describe the behavior and appearance of the chart. Since the Safety Outlier Explorer is a Webcharts chart object, many default Webcharts settings are set in the <u>defaultSettings.js file</u> as <u>described below</u>. Refer to the <u>Webcharts documentation</u> for more details on these settings.

In addition to the standard Webcharts settings several custom settings not available in the base Webcharts library have been added to the Safety Outlier Explorer to facilitate data mapping and other custom functionality. These custom settings are described in detail below. All defaults can be overwritten by users.

Renderer-specific settings

The sections below describe each safety-outlier-explorer setting as of version 2.3.0.

settings.id_col

string

unique identifier variable name

default: "USUBJID"

settings.time_cols

array

visit metadata

```
{
    "type": "ordinal",
    "value_col": "VISIT",
    "label": "Visit",
    "order_col": "VISITNUM",
    "order": null,
    "rotate_tick_labels": true,
    "vertical_space": 100
 },
  {
    "type": "linear",
    "value_col": "DY",
    "label": "Study Day",
    "order_col": "DY",
    "order": null,
    "rotate_tick_labels": false,
    "vertical_space": 0
 }
settings.time_cols[].type
string
Visit variable data type
default: "ordinal"
settings.timecols[].valuecol
string
Visit variable name
```

default:

```
default: "VISIT"
settings.time_cols[].label
string
Visit variable label
default: "Visit"
settings.timecols[].ordercol
string
Visit ordering variable name
default: "VISITNUM"
settings.time_cols[].order
array
Visit order
default: none
settings.timecols[].rotatetick_labels
boolean
Rotate tick labels 45 degrees?
default: true
```

settings.time*cols[].vertical*space

number

Rotated tick label spacing

default: 100

settings.measure_col

string

measure variable name

default: "TEST"

settings.unit_col

string

measure unit variable name

default: "STRESU"

settings.value_col

string

result variable name

default: "STRESN"

settings.normal*col*low

string

LLN variable name

default: "STNRLO"

settings.normal*col*high

string

ULN variable name

default: "STNRHI"

settings.start_value

string

value of measure to display initially

default: none

settings.filters

array

an array of filter variables and associated metadata

default: none

settings.filters[].value_col

string

Variable name

default: none

settings.filters[].label

string

Variable label

default: none

settings.details

```
array
an array of ID-level variables and associated metadata
default:
{
    "value_col": "AGE",
    "label": "Age"
 },
  {
    "value_col": "SEX",
    "label": "Sex"
 },
    "value_col": "RACE",
    "label": "Race"
 }
settings.details[].value_col
string
Variable name
default: "AGE"
settings.details[].label
```

string

```
Variable label
default: "Age"
settings.tooltip cols
object
added in v2.4. Takes an array of col/label objects and adds those values to the
end of the standard point tooltips. The following would add Date to the tooltip:
{tooltip_cols: [{label:"Date",value_col:"DT"}]};
default: null
settings.multiples_sizing
object
width and height of small multiples
settings.multiples_sizing.width
number
Width
default: 300
settings.multiples_sizing.height
number
Height
```

default: 100

settings.normalrangemethod

string

method for identifying the normal range

default: "LLN-ULN"

settings.normal*range*sd

number

the number of standard deviations from the mean with which to define the normal range

default: 1.96

settings.normal*range*quantile_low

number

this quantile defines the lower bound of the normal range

default: 0.05

settings.normal*range*quantile_high

number

this quantile defines the upper bound of the normal range

default: 0.95

settings.visitswithoutdata

boolean

controls display of visits without data for the current measure

default: false

settings.unscheduled_visits

boolean

controls display of unscheduled visits

default: false

settings.unscheduled*visit*pattern

string

a regular expression that identifies unscheduled visits

default: "/unscheduledlearly termination/i"

settings.unscheduledvisitsvalues

array

an array of strings that identify unscheduled visits; overrides unscheduled visit pattern

default: none

settings.line_attributes

object

an object that defines the line color, thickness, and opacity

settings.line_attributes.stroke

string

Line Color

default: "black"

settings.line_attributes.stroke-width

number

Line Thickness

default: 0.5

settings.line_attributes.stroke-opacity

number

Line Opacity

default: 0.75

settings.point_attributes

object

an object that defines the point color, radius, and opacity and its outline color, thickness, and opacity

settings.point_attributes.stroke

string

Point Outline Color

```
default: "rgb(102,194,165)"
settings.point_attributes.stroke-width
number
Point Outline Thickness
default: 0.5
settings.point_attributes.stroke-opacity
number
Point Outline Opacity
default: 1
settings.point_attributes.fill
string
Point Color
default: "rgb(102,194,165)"
settings.point_attributes.radius
number
Point Radius
default: 3
settings.point_attributes.fill-opacity
```

number

Point Opacity

default: 1

{

Webcharts settings

The object below contains each Webcharts setting as of version 2.3.0.

```
x: {
    column: null, //set in syncSettings()
    type: null, //set in syncSettings()
    behavior: 'raw'
},
y: {
    column: null, //set in syncSettings()
    stat: 'mean',
    type: 'linear',
    label: 'Value',
    behavior: 'raw',
    format: '0.2f'
},
marks: Γ
    {
        per: null, //set in syncSettings()
        type: 'line',
        attributes: {
             'clip-path': 'url(#1)'
        },
        tooltip: null //set in syncSettings()
    },
    {
        per: null, //set in syncSettings()
        type: 'circle',
        attributes: {
```

```
'clip-path': 'url(#1)'
     },
     tooltip: null //set in syncSettings()
     }
],
resizable: true,
margin: { right: 20 }, //create space for box plot
aspect: 3
```

Testing Logs

Interactive log of QC for all code updates. Includes code reviews, feature testing and regression testing for all releases. see original

PullRequest	Requester	DateRequested		Branch
0007	jwildfire	2016-04- 15T16:34:39Z	adjust-ticks	2. 4.16.1
0014	jwildfire	2016-05- 25T19:55:47Z	reorg	
0015	nbryant	2016-06- 14T17:17:40Z	build-fix	
0026	jwildfire	2016-12- 17T16:15:37Z	v1.3.0	
0028	samussiah	2017-01- 09T17:39:14Z	custom-filters	
0028	samussiah	2017-01- 09T17:39:14Z	custom-filters	

0028	samussiah	2017-01- 09T17:39:14Z	custom-filters
0029	samussiah	2017-01- 09T20:06:04Z	participant-details
0030	samussiah	2017-01- 10T15:39:40Z	small-multiple-layout
0031	samussiah	2017-01- 10T16:53:59Z	consistent-axes
0031	samussiah	2017-01- 10T16:53:59Z	consistent-axes
0031	samussiah	2017-01- 10T16:53:59Z	consistent-axes
0031	samussiah	2017-01- 10T16:53:59Z	consistent-axes
		2017-01-	

0032	samussiah	2017-01- 10T19:20:42Z	rotate-x-tick-labels
0032	samussiah	2017-01- 10T19:20:42Z	rotate-x-tick-labels
0033	samussiah	2017-01- 10T21:31:52Z	sample-population
0033	samussiah	2017-01- 10T21:31:52Z	sample-population
0033	samussiah	2017-01- 10T21:31:52Z	sample-population
0033	samussiah	2017-01- 10T21:31:52Z	sample-population
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks

10T19:20:42Z rotate-x-tick-labels

0032

samussiah

0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0034	samussiah	2017-01- 17T20:19:05Z	custom-marks
0038	samussiah	2017-01- 20T22:26:39Z	small-multiple-fixes
0038	samussiah	2017-01- 20T22:26:39Z	small-multiple-fixes
0044	samussiah	2017-02- 02T15:26:52Z	y-axis-label-fix
0045	jwildfire	2017-02- 23T21:46:07Z	master
0047	jwildfire	2017-03- 07T02:34:13Z	v2.0.1
0052	jwildfire	2017-10- 12T22:53:20Z	v2.0.2-dev
0055	jwildfire	2017-10- 16T20:12:56Z	v2.0.3-dev
0059	jwildfire	2017-11- 17T03:03:17Z	y-axis-controls

0059	jwildfire	2017-11- 17T03:03:17Z	y-axis-controls
0060	jwildfire	2017-11- 17T03:04:07Z	v2.1.0-dev
0061	jwildfire	2017-11- 17T03:49:06Z	check-values
0065	samussiah	2017-12- 05T18:00:02Z	y-axis-controls
0073	samussiah	2018-01- 26T20:03:40Z	v2.2.0-dev
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01-	visit-order

		26T20:04:24Z	
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
0074	samussiah	2018-01- 26T20:04:24Z	visit-order
		20120.01.212	
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
		2019 01	
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits
0076	samussiah	2018-01- 26T21:48:00Z	unscheduled-visits

0077	samussiah	2018-01- 26T22:04:27Z	visits-without-data
0077	samussiah	2018-01- 26T22:04:27Z	visits-without-data
0078	samussiah	2018-01- 26T22:27:45Z	settings-schema
0078	samussiah	2018-01- 26T22:27:45Z	settings-schema
0081	samussiah	2018-02- 01T22:14:05Z	revert-78-settings-schema
0082	samussiah	2018-02- 01T22:15:02Z	revert-77-visits-without-data
0083	samussiah	2018-02- 01T22:18:35Z	visits-without-data
0084	samussiah	2018-02- 01T22:20:20Z	settings-schema
0085	samussiah	2018-02- 02T22:07:32Z	v2.2.1-dev
0087	samussiah	2018-02- 27T14:56:09Z	v2.2.3-dev
0087	samussiah	2018-02- 27T14:56:09Z	v2.2.3-dev
0089	samussiah	2018-03- 27T13:39:27Z	v2.3.0-dev
0092	samussiah	2018-03- 28T15:48:25Z	darker-lines
0092	samussiah	2018-03- 28T15:48:25Z	darker-lines

0092	samussiah	2018-03- 28T15:48:25Z	darker-lines
0093	samussiah	2018-03- 30T18:51:27Z	normal-range-controls
0093	samussiah	2018-03- 30T18:51:27Z	normal-range-controls
0094	samussiah	2018-03- 30T19:24:53Z	fix-tick-rotation
0094	samussiah	2018-03- 30T19:24:53Z	fix-tick-rotation
0095	samussiah	2018-03- 30T20:39:52Z	append-units
0095	samussiah	2018-03- 30T20:39:52Z	append-units
0098	samussiah	2018-04- 02T15:15:40Z	add-test-page-script
0102	brittsikora	2018-04- 25T20:47:33Z	brittsikora-patch-1
0104	rtbailey	2018-08- 09T19:44:26Z	rtbailey-patch-1
0108	jwildfire	2018-10- 03T20:51:30Z	v2.4.0-dev
0108	jwildfire	2018-10- 03T20:51:30Z	v2.4.0-dev
0108	jwildfire	2018-10- 03T20:51:30Z	v2.4.0-dev
0110	jwildfire	2018-10- 04T04:20:15Z	sendToBackOnClick
0110	jwildfire	2018-10- 04T04:20:15Z	sendToBackOnClick
0110	jwildfire	2018-10- 04T04:20:15Z	sendToBackOnClick
0112	samussiah	2018-10- 10T21:06:13Z	fix-mark-rotation-to-front

0112	samussiah	2018-10- 10T21:06:13Z	fix-mark-rotation-to-front
0112	samussiah	2018-10- 10T21:06:13Z	fix-mark-rotation-to-front
0112	samussiah	2018-10- 10T21:06:13Z	fix-mark-rotation-to-front
0112	samussiah	2018-10- 10T21:06:13Z	fix-mark-rotation-to-front
0115	pburnsdata	2019-01- 08T21:14:35Z	pburnsdata-patch-1
0116	samussiah	2019-02- 07T15:33:48Z	dev-v2.5.0
0117	samussiah	2019-02- 07T16:56:51Z	invalid-starting-measure
0117	samussiah	2019-02- 07T16:56:51Z	invalid-starting-measure
0118	samussiah	2019-02- 07T18:35:28Z	extend-y-domain
0118	samussiah	2019-02- 07T18:35:28Z	extend-y-domain
0119	samussiah	2019-02- 07T23:06:45Z	add-y-increment
0119	samussiah	2019-02- 07T23:06:45Z	add-y-increment
0120	samussiah	2019-02- 07T23:09:39Z	removed-records
0120	samussiah	2019-02- 07T23:09:39Z	removed-records
0121	samussiah	2019-02- 12T21:35:52Z	custom-mark-highlighting
0121	samussiah	2019-02-	custom-mark-highlighting

		12T21:35:52Z	
0121	samussiah	2019-02- 12T21:35:52Z	custom-mark-highlighting
0123	jwildfire	2019-03- 25T18:48:02Z	dev-v2.5.1
0124	jwildfire	2019-03- 27T19:40:42Z	dev-v2.5.2
0126	samussiah	2019-04- 09T13:53:10Z	dev-v2.5.3
0133	jwildfire	2019-05- 15T23:32:11Z	dev-v2.5.4
0133	jwildfire	2019-05- 15T23:32:11Z	dev-v2.5.4
0133	jwildfire	2019-05- 15T23:32:11Z	dev-v2.5.4
0134	dependabot[bot]	2019-07- 12T19:32:52Z	dependabot/npm_and_yarn/lo 4.17.14
0137	jwildfire	2019-08- 09T19:51:10Z	dev-v2.5.5
0138	jwildfire	2019-08- 09T19:51:44Z	addCustomEvent
0139	jwildfire	2019-08- 14T05:57:18Z	dev-v2.6.0
0140	jwildfire	2019-08- 14T05:57:59Z	noverlap
0142	jwildfire	2019-08- 19T22:10:24Z	hoverlap
0144	samussiah	2020-01- 31T16:28:15Z	dev-v2.6.1